



CE 2777

USER INFORMATION

These products are classed as Cat.III Personal Protective Equipment (PPE) by the European PPE Regulation (EU) 2016/425 and have been shown to comply with this Regulation through the Harmonized European Standard(s) EN ISO 374-1:2016+A1:2018, EN ISO 21420:2020, EN ISO 374-5:2016.

Manufacturer: Pronto Direct Limited
Address: 8 Hornsby Square, Southfields Business Park, Basildon, Essex SS15 6SD, UK

Product reference:	Description:
FL4001-FL4005	Disposable Powder-Free Nitrile Gloves Blue(XS-XL)
FL4051-FL4055	Violet Blue(XS-XL)
FL4041-FL4045	Cobalt Blue (XS-XL)
FL4021-FL4025	White(XS-XL)
FL4011-FL4015	Black(XS-XL)
FL4032-FL4034	Pink (S-L)
FL4072-FL4074	Long Cuff White(S-L)
FL4062-FL4065	Long Cuff Blue(S-XL)
FL4082-FL4085	Long Cuff Violet Blue(S-XL)

Size: 1-5 (XS - XL)

Notified Body responsible for certification and ongoing conformity:

SATRA Technology Europe Ltd,
Bracetown Business Park,
Clonee, Dublin 15, D15 YN2P Ireland
Notified Body number: 2777



LOT



Classification:

EN ISO 374-1:2016+A1:2018 /Type B	Level	EN ISO 374-4:2013 Degradation%
40% Sodium Hydroxide (K)	6	1.4
30% Hydrogen Peroxide (P)	6	24.0
37% Formaldehyde (T)	6	10.9

EN ISO 374-1/Type B



KPT

EN ISO 374-5:2016	Level
Protection against Bacteria and Fungi	PASS

EN ISO 374-5:2016

Protection against Viruses

PASS



Virus

Note: The penetration resistance has been assessed under laboratory conditions and related only to test specimen.

Permeation levels are based on breakthrough times as follows:

Permeation performance level	1	2	3	4	5	6
Measured breakthrough time (min)	>10	>30	>60	>120	>240	>480

1. Intended Use:

The protective gloves (Disposable Nitrile Gloves, Powder Free) against chemical risks and microorganism are used for general purpose, food service, cleaning room and similar use. They are worn on the users' hands to protect their hands and fingers against chemical harm and microorganism contamination.

All materials used in the product: WILCAST, Willine-123-3, WILBELOW, Calcium Carbonate, Anti-webbing agent, Hydroxyethyl Cellulose, Alcohols C16-18 Ethoxylated.

The glove does not contain any substances that are known to cause allergies.

The gloves are for single use only.



The labels are printed in each box

The Gloves have no mechanical protection offered.

2. Performance and limitation of use:

This product has been tested in accordance with EN ISO 374-1:2016+A1:2018, EN ISO 21420:2020 and EN ISO 374-5:2016 and achieved the following performance levels.

2.1 EN ISO 374-1:2016+A1:2018

- This information does not reflect the actual duration of protection in the workplace and the difference between mixtures and pure chemicals.
- The chemicals resistance has been assessed under the laboratory conditions from samples taken from the palm only (except in cases where the gloves is equal to or over 400mm-where the cuff is tested also) and relates only to the chemicals tested. It can be different if the chemical is used in a mixture.
- It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation.
- When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistance gloves.
- Before usage, inspect the gloves for any defect or imperfections, do not use damaged gloves.

2.2 EN ISO 21420:2020

Size: the sizes of the gloves are XS, S, M, L. The equivalence to the sizes given in EN ISO 21420:2020 are shown in the below table:

Glove size	EN ISO 21420:2020
XS	6
S	7
M	8
L	9
XL	10

2.3 EN ISO 374-2:2019

This product has been tested for air leaks and water leaks and have been found to meet the pass level required by the standard EN ISO 374-2:2019.

2.4 EN ISO 374-4:2019

Degradation results indicate the change in puncture resistance of the gloves after exposure to the challenge chemical.

3. Protection limited:

No gloves may be used as protection for all chemicals. A glove may protect against a specific chemical, but it may not protect the wearer from another. If a glove protects the wearer, it will not protect the wearer forever, as the glove material will deteriorate. Therefore, the following must be considered when choose which gloves to be worn to protect against chemical exposure.

4. Other information:

4.1 Factors to consider when choose gloves:

Chemical to be used: consult the compatibility charts to ensure that the gloves will protect you.

Dexterity needed: the thicker the glove, typically the better the chemical protection.

Extent of the protection required: determine if a wrist length glove provides adequate protection.

4.2 Rules for glove use in the medical field:

Wear the correct gloves when needed.

Wear gloves no longer than 2 hours.

Wash hands once globes have been removed.

Disposable gloves must be discarded once removed. Do not save for future use.

Remove gloves before touching personal items, such as phones, computers, pens and one's skin.

Do not wear gloves out of working office. If gloves are needed to transport anything, wear one glove to handle to transported item.

If for any reason a glove fails, and chemicals come into contact with skin, consider it an exposure and seek medical attention.

4.3 Find the right size:

Find the appropriate size of the gloves. User can determine if their gloves too small by extending their fingers until straight.

4.4 Donning:

1. Remove all hand and wrist jewelry, and wash the hands before donning.
2. Place the gloves on the prepared work surface.
3. The user puts a glove on his/her dominant hand by grabbing it with the other hand, remembering to only touch the inside of the gloves, and slipping it over the dominant hand until it reaches finger level.
4. The wearer uses the gloved dominant hand to slip the other glove onto the non-dominant hand.
5. Once both glove are on, the users can touch the outside of the gloves to ensure a proper fit.

4.5 Doffing:

1. Using the dominant hand, users start by grabbing the outside of the glove on the non-dominant hand on the palm side near the cuff.
2. Pull the glove off the non-dominant hand and place it in the gloved hand, balling it up.
3. Slip two fingers under the cuff of the other hand glove and carefully peel it off the hand without touching the wrist, turning the remaining glove inside out as it is removed and in turn encasing the first glove.
4. The gloves can be disposed.

4.6 Storage and transport:

When not in use, store the products in a well-ventilated area away from extremes of temperature.

4.7 Packaging

100pieces/Box & 10 Box/Carton

4.8 Obsolescence

These gloves have a limited lifetime and should not be kept more than 5 years before use.

4.9 EU Declaration of Conformity

The EU Declaration of Conformity will follow the product to the market.